




# PORTABLE ACCESS RAMP SYSTEMS



**martyn industrials**  
**access ramp solutions**  
[www.martyn-industrials.com](http://www.martyn-industrials.com)



# A unique aid to personal mobility!

A simple but sophisticated folding ramp design, manufactured using high-tech materials which results in an easy to handle ramp that can be deployed **and** used by most people in seconds, with a minimum of experience or training.



## DEC PAC® Senior

The **Invacare Senior** ramp is designed without edge barriers, to be used where the doorway is narrower than normal (less than 800mm), or to offer a less expensive ramp solution for loading unoccupied powered wheelchairs into the rear of a Car or Van.

## DEC PAC® EBL

The **Invacare EBL** (Edge Barrier Limiting) model has been developed and tested with the assistance of Wheelchair users and healthcare professionals to give wheelchair users more Independent access to shops, offices, government buildings and public transport.

The Edge Barrier (50mm high) is designed to prevent wheelchairs steering off the edge of the ramp, especially important on the longer models.

## DEC PAC® Personal

The **Personal** is a small ramp which can be stored in a carrying bag.

The low weight enables it to be carried on your shoulder (like a bag) or hooked on the back of a rucksack.

This portability means it can be used wherever easy access is required, at home, at the workplace, on the street to overcome an entrance step!



**Invacare 120 Senior** and **120 EBL** ramps are now available in WIDE (86cm), to meet current EU DDA requirements for ramps used in transport applications.



by



**martyn industrials**  
access ramp solutions



## **DECPAC®** **Multi-Purpose**

Small lightweight  
supplied with its own

can be carried  
(as a rucksack) or  
attached to a wheelchair!

It can be used  
anywhere it is required - in the  
office or even in the  
road kerb or shop

The **Multi-Purpose** ramp is a longer, wider variation of the **Personal** model, but its light weight makes it just as portable.

The extra length and width make this model more versatile and can enable wheelchair users to climb larger steps and higher kerbs.

## **DECPAC®** **Carrying bag**

This purpose-made carrybag is designed to carry **Personal** and **Multi-Purpose** ramps securely and safely, either on your shoulder or attached to the wheelchair.

The heavyweight waterproof material helps to ensure the ramp stays dry and in good condition for when it is required.

## RAMP MODEL SPECIFICATION

Model	Width <i>largeur/breit</i>	Length <i>longueur/lange</i>	Weight <i>poids/gewicht</i>	Max. Step Height <i>hauteur max'm/ max.höhe</i>
-------	-------------------------------	---------------------------------	--------------------------------	--

### WHEELCHAIR RAMPS

<b>Personal</b>	74 cm	70 cm	4 kg	17.5 cm
<b>Multi-Purpose</b>	83 cm	87 cm	6 kg	22 cm

### PORTABLE RAMPS

Model	Senior	EBL	All Models	Senior/Standard	EBL	
<b>900</b>	75 cm	82 cm	90 cm	7 kg	8 kg	24 cm
<b>120</b>	75 cm	82 cm	120 cm	8 kg	9 kg	30 cm
<b>135</b>	75 cm	82 cm	135 cm	8.5 kg	11 kg	34 cm
<b>150</b>	75 cm	82 cm	150 cm	11 kg	12 kg	37.5 cm
<b>165</b>	75 cm	82 cm	165 cm	12 kg	13 kg	41 cm
<b>200</b>	75 cm	82 cm	200 cm	13 kg	16 kg	50 cm

### TRANSPORT RAMPS DDA COMPLIANT CE/TRANS/505/REV2/R107

<b>900 Wide</b>	86 cm	90 cm	N/A	8.5 kg	24 cm
<b>120 Wide</b>	86 cm	120 cm	8.5 kg	9.5 kg	30 cm

**ALL ramps have a Safe Working Load (SWL) of 300 kg**



### REFURBISHMENT AND SUPPORT

Ramps can be refurbished and/or repaired in the event of wear and damage. Please contact us for details.

### AS SUPPLIED TO



**martyn industrials**  
access ramp solutions

**For more information please contact:**

Gordon Keen on **07850 856201** or Mike Houston on **07850 113498**

5 Brunel Way, Durranshill Industrial Estate, Harraby, Carlisle CA1 3NQ

📞 01228 544 000 📠 01228 544 001 📧 enquiries@martyn-industrials.co.uk

**www.martyn-industrials.com**